Names:	Section_	
Chem 226/ Fall 2009		Dr. Rusay
	Optical Rotation / Polarimetry (PART 1)	

Refer to Lab Text/Guide pp.54-57, Class Notes, and the Web page listed below

Work in pairs. There will be no need for a formal report and records in your lab notebook for this experiment. There are 3 PARTS. Complete this form and the questions in PART 2 for the on-line *jmol* structures, and turnin before beginning PART 3 (Analysis of the Experimental Resolution of (+/-) Ibuprofen).

PART I: Optical rotation, optical purity, enantiomeric excess

http://chemconnections.org/organic/chem226/226assign-09.html#polarimetry

Olfactory discrimination of enantiomers is possible as you have seen with carvone. In PART I of this experiment you will determine the optical purity of a sample of carvone using a polarimeter and relate your experimental results to the enantiomer's smell, physical properties and absolute configuration. Each partner should independently determine α for the unknown carvone solution A that has been prepared for you. Unknown B has it's data provided. Take the average of the two and then calculate $[\alpha]$ for each of the carvone unknowns. Show your calculations below the Table.

Experimental Data:

Cell path length =	Temperature =	$\lambda = 589 nm$	solvent =	$\alpha_{\text{solvent}} =$
100. mm	25 °C	(sodium D)	ethanol	0°

	Mass	Volume	α_1	α_2	$\alpha_{ m avg}$	[α]	Smell: (mint,
	(mg)	(mL)				(calc.)	caraway or
							cannot tell)
Unknown A	4,002	25.00					
Unknown B	3,945	25.00	+7.0°	+7.1°	+7.05°		

Calculations:

$$\left[\alpha\right]_{A} = \left[\alpha\right]_{B} =$$

Consult the chemical literature and complete the following table of physical / optical data for the respective carvone enantiomers. (To determine the absolute configurations (R- or S-) refer to the structure below of d-carvone.)

	boiling point	density	[α]	Abs. Config.
d-carvone				
l-carvone				

Using the literature and experimental data complete the following questions for the unknowns A and B. Show your calculations for optical purity and enamtiomeric excess.

	optical purity	% R-	% S-	Enantiomeric Excess: (%) indicate d- or l-	Smell
Unknown					
I I I I I I I I I I I I I I I I I I I					
Unknown B					

	Unknown					
	A Unknown					
	В					
Calculat Optical P	tions: turity (Enantion	neric Exces	s) A			
Optical P	turity B (Enanti	omeric Exc	ess B)			
A) Are th	ne observed sma	ells consiste	ent with you	ur results? B	riefly discuss.	
					g . 1) Is carvone co Briefly explain you	2) Would