


Chem 108: Lab Week 13

Sign in / Pick up Papers

Pick a vial and a plastic
dropper

Using the vial number, sign-
in next to your name on
the Lab roster



Chem 108: Class/ Lab Week 13

TODAY:

1) Fluid Exchange
(Handout) *Due Next Lab*

2) To Do: Acid-Base Equilibrium Experiment
(Handout)

Data completed & signed before leaving Lab

3) To Do: Unknown acid titration

Data completed & signed before leaving Lab

Follow Instructions

<http://chemconnections.org/general/chem120/fluid-ex.108.html>


Chem 108: Class/ Lab Week 13

TODAY:

Fluid Exchange (Handout)


3) You have been assigned a geographical location for your Global Residence. Check the *Global Homelands Map*, which follows, for your location and if necessary move to your place of residence.

<http://chemconnections.org/general/chem120/fluid-ex.108.html>



Wait for Dr. R's instructions on exchanging fluids,
keeping records, and using the forms provided.

Acid-Base Indicators



Indicators

Natural Indicators



Acid-Base Indicators

	pH range for color change													
	0	2	4	6	8	10	12	14						
Methyl violet	Yellow							Violet						
Thymol blue		Red				Yellow			Yellow			Blue		
Methyl orange			Red											
Methyl red				Red										
Bromthymol blue					Yellow							Blue		
Phenolphthalein							Colorless					Pink		
Alizarin yellow R											Yellow			Red

Chem 108: Class/ Lab

Week 13

- 1) Fluid Exchange
(Handout) *Due Next Lab*

TODAY:

- 2) To Do: Acid-Base Equilibrium Experiment
(Handout)

Data table completed & signed before leaving Lab


Fully Completed Handout Due Next Week

Lab: universal Indicator, litmus & pH paper



	Red Litmus	Blue Litmus	Solution pH	Description
A	red	red	1	acid
B	blue	blue	13	base
C	red	blue	7	neutral
D	red	blue	7	neutral
E	blue	blue	10	base
F	blue	blue	10	base
G	blue	blue	10	base
H	blue	blue	10	base
I	blue	blue	10	base


J	NaCl(aq)	red	red	7.0
K	NaOH(aq)	blue	blue	10.9
L	Mg(OH) ₂	blue	blue	12.2
M	Orange juice	red	red	3.8
N	Milk	red	red	6.4
O	Saliva (spit) and blood	blue	blue	7.3
P	Vinegar	red	red	1.9
Q	Butter (pH 7)	red	blue	7.0



	Red Litmus	Blue Litmus	Solution pH	Description
A	red	red	1	acid
B	blue	blue	13	base
C	red	blue	7	neutral
D	red	blue	7	neutral
E	blue	blue	10	base
F	blue	blue	10	base
G	blue	blue	10	base
H	blue	blue	10	base
I	blue	blue	10	base

J	NaCl(aq)	red	red	7.0
K	NaOH(aq)	blue	blue	10.9
L	Mg(OH) ₂	blue	blue	12.2
M	Orange juice	red	red	3.8
N	Milk	red	red	6.4
O	Saliva (spit) and blood	blue	blue	7.3
P	Vinegar	red	red	1.9
Q	Butter (pH 7)	red	blue	7.0

Lab pH: pH Meter (Chem 120)



	Red Litmus	Blue Litmus	Solution pH	Description
A	red	red	1.0	acid
B	blue	blue	13.0	base
C	red	blue	7.0	neutral
D	red	blue	7.0	neutral
E	blue	blue	10.1	base
F	blue	blue	10.1	base
G	blue	blue	10.1	base
H	blue	blue	10.1	base
I	blue	blue	10.1	base

J	NaCl(aq)	red	red	7.0
K	NaOH(aq)	blue	blue	10.9
L	Mg(OH) ₂	blue	blue	12.2
M	Orange juice	red	red	3.8
N	Milk	red	red	6.4
O	Saliva (spit) and blood	blue	blue	7.3
P	Vinegar	red	red	1.9
Q	Butter (pH 7)	red	blue	7.0


<https://phet.colorado.edu/en/simulation/ph-scale-basics>

The pH Scale

pH: the negative logarithm
of the
hydrogen ion concentration.



Beer's pH ~ 4



Quantitative,
logarithmic, numeric
scale based on testing
the electric current of
aqueous solutions &
relating it to the
equilibrium
concentration of the
hydrogen ion,
[H⁺]_(aq) = [H₃O⁺]_(aq)

Introduced in 1909 by
Søren Sørensen, Danish
brewer/chemist, as a
convenient way of
expressing acidity.....
Providing much improved
quality control in brewing.

<http://www.chemconnections.org/general/chem108/Acids-Bases%20Guide.html>

Chem 108: Class/ Lab

Week 13

- 1) Fluid Exchange
(Handout) *Due Next Lab*
2) To Do: Acid-Base Equilibrium Experiment
(Handout)

Data completed & signed before leaving Lab

TODAY:

- 3) To Do: Unknown acid titration
Data completed & signed before leaving Lab

Chem 108: Lab

Standardization (Part 1) will not be done. Molarity of sodium hydroxide standard solution will be provided.

Name: _____
Section: _____

Individual Titrations of Unknown Acid To Do Today

Report Form - Acid Base Titration

Part 1-Standardization of NaOH Solution


Molarity of HCl used						
Titration	1	2	3	4	5	6
Base buret, final reading (mL)						
Base buret, initial reading (mL)						
Volume of base used (mL)*						
Molarity of NaOH (M)*						
Average molarity of NaOH*						M

Show calculations for each of the entries in the Data Table marked with * on the calculations page for one titration.

<https://www.youtube.com/watch?v=9Dk882xLvNE>

Acid-Base Titration

<https://www.youtube.com/watch?v=9Dk882xLvNE>



Part 2: Unknown Acid

Equipment

From the stockroom:
plastic 1 L bottle
50 mL buret
buret clamp
25 mL vol. pipet and bulb
From the common drawer:
ring stand
From your drawer:
funnel
125 mL flask
250 mL flask
2 beakers (one for waste)
wash bottle

Chem 108: Lab

Part 2:

To Do (individually) today. Have signed before leaving lab.

Part 2-Determination of Unknown Acid

Unknown code						
Average Molarity of Base from Part 1			0.2051 mol/L			
Titration	1	2	3	4	5	6
Base buret, final reading (mL)						
Base buret, initial reading (mL)						
Volume of base used (mL)*						
Molarity of unknown acid (M)*						
Average molarity of unknown (M)*						M

3 trials must be within +/- 0.20 mL



Show the calculations for each of the entries in the Data Table marked with * on the calculations page for one titration.

Take a clean, dry, 125 mL erlenmeyer flask to the stockroom window and get unknown acid solution. Record unknown number. Have data page signed before leaving lab today.

<https://www.youtube.com/watch?v=9Dk882xLvNE>


Neutralization Reactions

Titration



Neutralization Reactions



Titration



pH curves
Chem 120/121

Unkown Acid Neutralization

Net Ionic Equation/ Calculation




25.00 mL of $M_{H^+}(aq) = ?$ (unknown monoprotic nitric acid solution) was titrated with a sodium hydroxide solution, $M_{OH^-} = ? 0.2162 M$. It required 24.20 mL as an average of three trials which were within +/- 0.20 mL to reach a faint pink color.

$M_{H^+}(aq) = ?$

$$M_{H^+} = [M_{OH^-} \times V_{OH^-} / V_{H^+}] [? mol_{H^+} / ? mol_{OH^-}]$$

Unkown Acid Neutralization

Net Ionic Equation/ Calculation



25.00 mL of $M_{H^+ \text{ aq}} = ?$ (unknown monoprotic acid solution) was titrated with a sodium hydroxide solution, $M_{OH^-} = ? 0.2162 \text{ M}$. It required 24.20 mL as an average of three trials which were within $\pm 0.20 \text{ mL}$ to reach a faint pink color.


$$?M_{H^+} = [M_{OH^-} \times V_{OH^-} / V_{H^+}] [? \text{ mol}_{H^+} / ? \text{ mol}_{OH^-}]$$

$$= \frac{0.2162 \text{ mol}_{OH^-} \times 0.02420 \text{ L}_{OH^-} \times 1 \text{ mol}_{H^+}}{\text{L}_{OH^-} \times 0.02500 \text{ L}_{H^+} \times 1 \text{ mol}_{OH^-}} = 0.2093 \text{ M}_{H^+}$$

QUESTION

A 35.00 mL sample of a monoprotic acid of unknown concentration was titrated with 42.30 mL of 0.2250 M KOH. What is the concentration of the unknown acid?

- A. 0.0930 M
- B. 0.3030 M
- C. 0.2719 M
- D. 0.1356 M
- E. 0.3720 M




$$?M_{H^+} = [M_{OH^-} \times V_{OH^-} / V_{H^+}] [? \text{ mol}_{H^+} / ? \text{ mol}_{OH^-}]$$

QUESTION

A 35.00 mL sample of sulfuric acid (a di-protic acid) of unknown concentration was titrated with 42.30 mL of 0.2250 M KOH. What is the concentration of the unknown acid?

- A. 0.0930 M
- B. 0.3030 M
- C. 0.2719 M
- D. 0.1356 M
- E. 0.3720 M



$$? M_{H^+} = [M_{OH^-} \times V_{OH^-} / V_{H^+}] [? \text{ mol}_{H_2SO_4} / ? \text{ mol}_{OH^-}]$$

Chem 108: Lab Part 2: Week 13

To Do (individually) today

Part 2-Determination of Unknown Acid

Unknown code						
Average Molarity of Base from Part 1			0.2051 mol/L			
Titration	1	2	3	4	5	6
Base buret, final reading (mL)						
Base buret, initial reading (mL)						
Volume of base used (mL)*						
Molarity of unknown acid (M)*						
Average molarity of unknown (M)*					M	

Show the calculations for each of the entries in the Data Table marked with * on the calculations page for one titration.

Take a clean, dry, 125 mL erlenmeyer flask to the stockroom window and get unknown acid solution. Record unknown number. Have data page signed before leaving lab today.