

## Alcoholic Beverage Analysis

Purpose: The purpose is to determine the percent by volume of alcohol in an alcoholic beverage by distillation of a beer, ale, or wine, removing the water & ethanol.

- Procedure:
- ① Obtain beverage, twice amt. of what is to be used.
  - ② Beer = 50 mL, ale = 30 mL, table wine = 20 mL, & dessert wine = 10 mL.
  - ③ Shake until frothing stops if carbonated.
  - ④ Pipet approx. amt. of beverage into 500 mL dist. flask.
  - ⑤ Dilute 100 mL with d.i. water
  - ⑥ Add boiling chips, spray w/ anti-foam spray
  - ⑦ Assemble distillation apparatus.
  - ⑧ Heat gently, until frothing stops.
  - ⑨ Distill sample until 75 mL of distillate has been obtained in a 100 mL graduated cylinder.
  - ⑩ Transfer distillate to a 250 mL volumetric flask & dilute to volume with d.i. water.
  - ⑪ Mix well.
  - ⑫ Pipet 10 mL of contents of volumetric flask into 500 mL dist. fl.
  - ⑬ Pipet 10 mL of standard 0.2 M  $K_2Cr_2O_7$  into same flask.
  - ⑭ Record exact concentration of  $K_2Cr_2O_7$  solution.
  - ⑮ Add 50 mL of 3 M sulfuric acid to same flask.
  - ⑯ Add boiling chips & assemble reflux system.
  - ⑰ As it gets to a boil, adjust heat to a gentle boil, & continue refluxing for 90 minutes.
  - ⑱ While refluxing, do practice titrations.
  - ⑲ After 90 min turn off heat & cool.
  - ⑳ Transfer to 100 mL volumetric flask.
  - ㉑ Dilute & mix.
  - ㉒ Pipet 10 mL into 2 250 mL flasks, adding 5 drops of Ind., & titrate w/ 0.025 M iron(II) ammonium sulfate.
  - ㉓ Dissolve iron(II) ammonium sulfate with di  $H_2O$  in 250 mL fl.
  - ㉔ Calculations


Name

Course

Section

Lab partner

## Chemical Reactions:



## Data:

Alcohol: corona Extra Beer	
Alcohol Obtained:	50 mL
Mass $\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2$	2.4416 g
Titration #1 Reading	20.32 mL
Solution in Buret #1	50.00 mL
Actual Titrated Amount #1	29.68 mL
Titration #2 Reading	20.58 mL
Solution in Buret #2	50.00 mL
Actual Titrated Amount #2	29.42 mL
Titration #3 Reading	20.31 mL
Solution in Buret #3	50.00 mL
Titration #3 Reading Amount #3	29.67 mL
Average Titration	29.68 mL

## Calculations:

$$M_{\text{Fe}^{2+}} = \frac{2.4416 \text{ g Fe}^{2+}}{392.17 \text{ g Fe}(\text{NH}_4)_2(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}}$$

$$= 0.0249 \text{ mol/L Fe}^{2+}$$

$$M_{\text{Cr}_2\text{O}_7^{2-}} = \frac{0.0249 \text{ mol/L} \times 0.02968 \text{ L Fe}^{2+} \times \frac{1 \text{ mol Fe}^{2+}}{6 \text{ mol Cr}_2\text{O}_7^{2-}}}{0.01 \text{ L Cr}_2\text{O}_7^{2-}}$$

$$= 0.0123 \text{ mol/L Cr}_2\text{O}_7^{2-}$$

$$\text{moles Cr}_2\text{O}_7^{2-} (\text{excess}) = \frac{0.0123 \text{ mol Cr}_2\text{O}_7^{2-} \times 0.01 \text{ L Cr}_2\text{O}_7^{2-} \times 1000}{1 \text{ L}}$$

$$= 0.00123 \text{ moles Cr}_2\text{O}_7^{2-} (\text{excess})$$

$$\text{or } 1.23 \times 10^{-3} \text{ moles Cr}_2\text{O}_7^{2-} (\text{excess})$$

$$\text{moles Cr}_2\text{O}_7^{2-} (\text{total}) = 0.2 \text{ M Cr}_2\text{O}_7^{2-} \times 0.01 \text{ L Cr}_2\text{O}_7^{2-}$$

$$= 0.002 \text{ moles Cr}_2\text{O}_7^{2-}$$

$$\text{moles Cr}_2\text{O}_7^{2-} (\text{reacted w/ETOH}) = 0.002 \text{ moles Cr}_2\text{O}_7^{2-} - 0.00123 \text{ mol Cr}_2\text{O}_7^{2-}$$

$$= 7.7 \times 10^{-4} \text{ moles Cr}_2\text{O}_7^{2-} (\text{reacted})$$

$$\text{moles ETOH} (\text{reacted}) = 7.7 \times 10^{-4} \text{ mol Cr}_2\text{O}_7^{2-} \times \frac{3 \text{ mol ETOH}}{2 \text{ mol Cr}_2\text{O}_7^{2-}}$$

$$= 0.001155 \text{ moles ETOH} \times 25$$

$$= 0.0289 \text{ moles ETOH}$$

$$\text{grams ETOH} = 0.0289 \text{ mol ETOH} \times \frac{46.07 \text{ g ETOH}}{1 \text{ mol ETOH}}$$

$$= 1.33 \text{ g ETOH}$$

$$\text{WT:VOL} = \frac{1.33 \text{ g ETOH}}{50 \text{ mL}} \times 100 = 2.66\%$$

Name

Course

Section

Lab partner

Calculations:

$$\text{VOL:VOL}\% = \frac{1.33\text{g ETOH}}{(50\text{mL} \times 0.789)} \times 100\% = 3.37\%$$

$$\text{proof} = 3.37\% \times 2 \\ = 6.74 \text{ proof}$$

Conclusion:

I conducted 3 separate titrations, but I took the two that were within one percent of each other and took the average of those. It took 29.68mL for the solution to react. After conducting my calculations I was able to obtain an alcohol % of 3.37% and a proof of 6.74. The alcohol I used was regular Corona Extra Beer, which I originally obtained 50mL of. After researching the alcohol % of Corona Extra Beer, I discovered it has a 4.6% alcohol percent. The reason I believe why my calculated percent was lower was due to error. This experiment took 2 weeks to conduct & included lots of steps. The possibilities of error could be endless. Some prime examples would be improper titration, not refluxing long enough, & an incorrect reading of the buret. Although my final percent was not correct, I believe the tying together of ALL equations we have learned was the most vital part of this experiment.

Average of titrations:

$$\#1 \ 50.00\text{mL} - 20.32\text{mL} = 29.68\text{mL used}$$

$$\#3 \ 50.00\text{mL} - 20.33\text{mL} = 29.67\text{mL used}$$

$$29.68\text{mL} + 29.67\text{mL} = \frac{59.35\text{mL}}{2} = 29.675\text{mL}$$

$$\text{average titration} = 29.68\text{mL}$$

Grams of  $\text{Fe}(\text{NH}_4)_2(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O}$ :

$$392.17\text{g Fe}(\text{NH}_4)_2(\text{SO}_4)_2 \cdot 6\text{H}_2\text{O} \times 2.50\text{L} \times 0.025\frac{\text{mol}}{\text{L}} \\ = 2.451\text{g iron(II) ammonium sulfate}$$