# Chem 108: Lab Week 4

Sign in: Roster @ front of lab Pick up replacement page for today's experiment.

### Last Week's Experiment: Metric Measurement Completed individual forms pp. 12-15 due today by the end of lab.

		Name:	
		Section:	
	Report Form - Met	ric Measurement	
ength and Area			
Equipment Drawer			
	cm		m*
2. Large Test Tube			
	cm	m*	mm*
3. Crucible			
	cm	m*	mm*
4. Page		Hart British	
	cm	m*	mm*
Area of Page (Show you	calculations on the la	ast page of the Repo	rt Form.)
Length	Width	Area in cm <sup>2</sup>	Area in mm² & m²
cm	cm	cm <sup>2*</sup>	m2*
			mm²*

Show the calculations for each of the entries in the Data Table marked with * on the calculations page.	
---	--

1. Largest Test Tube			
	mL		L
2. Crucible			
	mL		L.
3. Smallest Test Tube			
	mL		L.
4. 250 mL beaker-graduated cyl	inder readings		1791 14 1 1 1 1 1
mL		mL	mL
Total		12/19/18	
	mL*		L*
5. 250 mL beaker-measured as	a cylinder		10:01 1 4 19
Height	Diameter		Radius*
cm		cm	cm
	Volume		cm3*

1. Crucible		1 1 1 1 1 1 1 1
	g	mg*
2. Crucible Lid		
	g	mg*
3. Crucible and Lid		
	g	mg*
4. Sum of Crucible and Lid		
Crucible		9
Lid		9
Sum		g
5. Equipment Slip		
	g	mg*

Include calculations. Staple together with partner's form and turn in before the end of lab.

Most legible report pages on top.

### Exp. 1 – Metric Measurement


Example of an acceptable set of student data, conversions, and calculations. DO NOT COPY. Use as a guide.

Volume 1. Largest Test Tube 0.06/0L+1-0.0001 L 6/10ml 1-0.1ml ml 2. Crucible 31,0mL++ 0.1mLmL 0.0310 L + - 0.0001 L 3. Smallest Test Tube 10.05 mL +1-0.01 mL 0,04005 L +1-0.00001 L 4. 250 mL beaker-graduated cylinder readings 97.3 ml +-0.1mlml 92.9 ml +-0.1ml 761 +-0.1 ml ml 266.3 ml +1-0.1ml ml+ 2663L +1-0,0001L 5. 250 mL beaker-measured as a cylinder Radius\* 3.5 cm ++ 0.1 cm cm 7.0 cm +1- 0.1 cm 96 cm + 1-01 cm 330, 9' (m 3+1-01/cm3

Show the calculations for each of the entries in the Data Table marked with \* on the calculations page

# Metric Measurement Conversions & Uncertainty

4.3 cm +/- 0.1cm



4.25 cm +/- 0.05cm

$$4.25 \text{ cm} + -0.05 \text{ cm} = 0.0425 \text{ m} + -0.0005 \text{ m}$$

$$4.2 \text{ cm}^2 + -0.1 \text{cm}^2 = 4.2 \text{ cm}^2 \text{ x 1m}/100 \text{cm x 1m}/100 \text{cm}$$
  
=  $0.00042 \text{m}^2 + -0.00001 \text{ m}^2$ 

	Adapted from Workshop Chemistr
Name(s):	
Worksheet: Units, Meas	surements, & Conversions
	ohxVNTg&list=PL8dPuuaLjXtPHzzYuWy6fYI dex=2 (11:23 min/sec)
1. How many significant figures are there	in the following numbers?
a) 42,000. L b) 0.4	9010 g
c) 0.00130 s d) 40	5,700,000 km
Complete the table. Provide ordinary de of unit. The first line has been complete	cimal form or scientific notation and the type d as an example for <b>mass</b> .
Ordinary Decimal Form	Scientific Notation
0.683 kg (mass)	$6.83 \times 10^{-1} \text{ kg}$
1365 mL (	mL
( )	1.034 × 10 <sup>1</sup> m
0.00350 μs (	μs
( )	1.75 × 10 <sup>-3</sup> cm <sup>3</sup>
1,605,000 nm (	nm
the following values using the specified	as read from a speedometer when traveling a

Worksheet (Handout): Due Today. Collaboration is encouraged. Turn in one with the names of all contributors.

http://chemconnections.org/general/chem108/Math %20%26%20Measurements-WKS.f18.pdf

### **CHEM 108**

Experiment 2: Measuring Density

## Comparison of Relative Properties FUNCTIONS



#### Definition of FUNCTION

- 1 : professional or official position : OCCUPATION His job combines the functions of a manager and a worker.
- 2 : the action for which a person or thing is specially fitted or used or for which a thing exists : PURPOSE
- : any of a group of related actions contributing to a larger action; especially: the normal and specific contribution of a bodily part to the economy of a living organism • The function of the heart is to pump blood through the body.
- 4 : an official or formal ceremony or social gathering They went to several functions during their college reunion weekend.
- 5 a: a mathematical correspondence that assigns exactly one element of one set to each element of the same or another set
  - **b**: a variable (such as a quality, trait, or measurement) that depends on and varies with another height is a *function* of age; *also*: **RESULT** illnesses that are a *function* of stress
- 6 : characteristic behavior of a chemical compound due to a particular reactive unit; also : FUNCTIONAL GROUP
- 7 : a computer subroutine; specifically : one that performs a calculation with variables (see <sup>2</sup>VARIABLE 1a) provided by a program and supplies the program with a single result

Density is a function of an object's mass and volume.

## Comparisons of Relative Properties FUNCTIONS

density: mass: volume

Functions & variables can be described interchangeably between:

**Density** 

1) Images (Symbols)



- 3) Sounds (Voice) Mass Volume
- 4) Numbers
- 5) Tables
- 6) Graphs
- 7) Formulas

Density is a function of an object's mass and volume.

### **Density**



http://www.density.com/what.htm

Density = Mass / Volume [g/mL or g/cm<sup>3</sup>; g/L]

mass (m)

$$d = \frac{m}{V} = \frac{156 \text{ g}}{20.0 \text{ cm}^3} = 7.80 \text{ g/cm}^3$$

volume (V)


density (d)


# Experiment 2 Measuring Density

- ➤ **Goal:** To measure diameter and height of a metal cylinder and calculate the respective volume.
  - To plot mass versus volume and determine slope of "trendline(s)".
  - Using 2 different methods, to measure mass and volume of a liquid, and to determine which method results in higher-precision (most decimal places)
  - To measure mass and volume of a solid using instruments of different precision, and determine which method results in *highest*—*precision*
- Work with same lab partner(s) as the Metric Measurement Experiment
  - Be sure to write yours and partners' names ON both REPORT FORMS DUE Next Week


### Chem 108 Lab: Week 4

### Volumes of regular shapes











### Reminder:

## Measuring with a 100 mL Graduated Cylinder:



Measurement using a volumetric pipet:



### meniscus

Bottom of the meniscus MUST be at calibration line for accurate measurement

### > Add a column to pg. 22 & label both of them

Each partner to measure with both Caliper and Ruler & compare





Method 3-Measuring Volume with Calipers	CALIPER	RULER
Mass, metal		
Height, metal		
Diameter, metal		
Radius, metal*		
Volume, metal*		
Density, metal*		
Densities of the Unknowns: Part A- Average Density of Unknown Liquid from Density of Liquid (Graduated Cylinder)	m the two methods	:
Part A- Average Density of Unknown Liquid from	m the two methods	
Part A- Average Density of Unknown Liquid from	m the two methods	:
Part A- Average Density of Unknown Liquid from Density of Liquid (Graduated Cylinder)  Density of Liquid (Volumetric Pipet)	m the two methods	:
Part A- Average Density of Unknown Liquid from		
Part A- Average Density of Unknown Liquid from Density of Liquid (Graduated Cylinder) Density of Liquid (Volumetric Pipet)  Average Density of Unknown Liquid*		
Part A- Average Density of Unknown Liquid from Density of Liquid (Graduated Cylinder)  Density of Liquid (Volumetric Pipet)  Average Density of Unknown Liquid*  Part B- Average Density of Unknown Metal from		
Part A- Average Density of Unknown Liquid from Density of Liquid (Graduated Cylinder)  Density of Liquid (Volumetric Pipet)  Average Density of Unknown Liquid*  Part B- Average Density of Unknown Metal from Density of Metal (Graduated Cylinder)		

### > Replace Part C, pg. 19 with handout

http://chemconnections.org/general/chem108/Measuring%20Density-graphing.2017.pdf



Measuring Density (REPLACEMENT for Part-C page 19)

Part C-Mass Versus Volume

Select either Data Set for Metal A or Data Set for Metal B below to plot. You will do one and your partner will do the other Data Set.

#### Metal A

Volume	Mass (g)
(cm3)	
7.89	17.22
6.80	18.11
7.92	21.21
9.75	22.25
8.17	23.19
9.84	25.44
10.1	26.36
11.4	28.29
11.8	28.73
44.7	20.60

#### Metal B


Volume	Mass (g)
(cm3)	
6.63	53.90
6.45	54.53
6.42	57.15
6.61	57.34
7.97	69.15
8.40	69.43
7.98	72.24
9.65	84.84
9.65	86.14
9.84	87.67

- Calculate the respective densities for the metals using an average of the Mass and Volume for each metal. (Report the values in the table on this page and attach it along with your graphs to the REPORT FORM pages and turn in.)
- 2. Graph the Mass (vertical axis) vs. Volume (horizontal axis) for your Data Set (graph paper is in Appendix C). Scale your graph to use as much of the graph paper as possible. Provide a Title, label the axes (include units), show each data point and draw the best possible straight line balancing the differences in the distance of the points on each side of the line using a ruler. Draw the line through the y-axis. Calculate the slope of the line, slope (m) = Δy / Δx, which equals the metal's density. (Report the values below.)
- The respective metals are either aluminum, density = 2.64 g/cm³ or copper, density = 8.94 g/cm³. Using these as accepted values, identify A and B, then calculate and record the calculated % error in the experimental densities using the two methods.

	Α	В
Metal identified		
Density (g/cm³) averaged		
Error (%) averaged		
Density (g/cm <sup>3</sup> ) graphed		
Error (%) graphed		




#### Method 2:




# **Density**Mass & Volume

Numbers Tables Graphs



#### Metal Mass & Volume 2270.0 2500.0 2000.0 1587 Mass (grams) 1500.0 Volume 1000.0 Volume 586.0 Mass (mL) 500.0 82.2 117.6 0.0 **Unknown Samples 1-6**




ables, & Graphs

Plot data using blank graph paper in back of Lab Manual pp. 119 – 121

Data & instructions are provided on handout

http:// chemconnections.org/general/ chem108/Measuring Densitygraphing.2017.pdf

- Plot all data for either
   Metal A (one partner) or
   Metal B (other partner):
   Can be on the same
   graph paper or separate.
   If on the same paper be
   sure to consult with
   partner on the scales to
   use.
- Complete the bottom table of handout and attach to the Report Form to turn in.



### Equation of a line: $\Delta y = m\Delta x + b$

$$y = y$$
 axis  $m = slope$   $x = x$  axis  $b = y$ -intercept

We're plotting: 
$$Mass = y axis$$
  $Volume = x axis$ 

➤ How are mass and volume related?

$$\frac{\Delta \ mass}{\Delta \ Volume} = density$$


We can rearrange this as: mass = density(Volume)If we compare to equation of a line:

$$mass = density(Volume) + 0$$
  
 $\Delta y = m \quad \Delta x \quad + b$ 

Now, what does the slope of our trendline represent?

### Plotting data

- A good plot should have:
  - An appropriate scale
  - Labeled axes (with units)
  - An appropriate title



Use a ruler to estimate best-fit line for data that appears to group together (It is NOT connect-thedots)

### Experiment 2 – Measuring Density

- "Extrapolate your lines through the y (vertical) axis"
  - This means draw your line all the way to the y-axis


### Error / Accuracy

 $Percent Error = \frac{Experimental \ value \ - True \ value}{True \ value} \times 100$ 

4 experimental sets of data

True = 1.0 g/mL Exp. = 1.0 g/mL

True = 1.0 g/mL Exp. = 1.1 g/mL



Value = Average

True =1.0 g/mL Exp. = 1.0 g/mL

True =1.0 g/mL Exp. = 1.2 g/mL

### **QUESTION**

Rank the relative accuracy of the three sets of data: a), b) and c). The accepted value is 8.08 mL.

Average	Average	Average
a)	b)	c)
8.38	8.38	7.99
Error = ?	Error = ?	Error = ?
Standard	Standard	Standard
deviation	deviation	deviation
a)	b)	c)
+/- 0.91	+/- 0.03	+/- 0.05

A) Accuracy: a > c > b

B) Accuracy: b > c > a

C) Accuracy: c > a = b

D) Accuracy: a = b > c

### **QUESTION**

Rank the relative accuracy of the three sets of data: a), b) and c). The accepted value is 8.08 mL.

Average	Average	Average
a)	b)	c)
8.38	8.38	7.99
Error = ?	Error = ?	Error = ?
Standard	Standard	Standard
deviation	deviation	deviation
a)	b)	c)
+/- 0.91	+/- 0.03	+/- 0.05

A) Accuracy: a > c > b

Answer:

C) Accuracy: c > a = b

B) Accuracy: b > c > a

D) Accuracy: a = b > c

### Experiment 2 – Measuring Density


- "Extrapolate your lines through the y (vertical) axis"
  - ➤ This means draw your line all the way to the y-axis

**Error / Accuracy & Precision** 

$$Percent Error = \frac{Experimental \ value \ - True \ value}{True \ value} \times 100$$

True = 1.0 g/mLExp. = 1.0 + /- 0.1 g/mL


True =1.0 g/mL Exp. = 1.1 +/- 0.1 g/mL



True = 1.0 g/mL Exp. = 1.0 + - 0.2 g/mL

True =1.0 g/mL Exp. = 1.2 +/- 0.2 g/mL

### Different Volumetric Tools



A and D have the same level of precision.

True

The following is the correct order of increasing precision: D < C < B. False

### **QUESTION**

Rank the relative precision of the three sets of data: a), b) and c). The accepted value is 8.08 mL.

Average	Average	Average
a)	b)	c)
8.38	8.38	7.99
Standard	Standard	Standard
deviation	deviation	deviation
a)	b)	c)
±/- 0.91	<b>⊥/-</b> 0 03	<b>4/-005</b>


Answer:


A) Precision: a > c > b B) Precision: b > c > a

C) Precision: a = b > c D) Precision: a > b > c

### > Complete both columns pg. 22

### Each partner to measure with both Caliper and Ruler & compare






#### Using a Caliper

Come up with partner having names on both lab forms.

Dr. R will provide unknown liquid & unknown metals, caliper & ruler.



Measuring the cylinder





Read instructions carefully & discuss with partner before beginning.


Return empty liquid vials & metal cylinders when finished.


Liquid is salt solution that can be poured down drain.




### Experiment 2 – Measuring Density

## Complete and record all measurements pp. 20-22 today. Omit pg. 23 and Trial 2 pg.20







Have Dr. R. sign **individual** forms before leaving lab today.

Only your measurement data is due to be signed today.

Be certain of calculations, graphs & questions before leaving lab.

### For Next Week

- Check Calendar for assignments
- Complete density calculations, graphs & Report Form pp.20-22, & pp. 24-25; attach completed replacement pg. 19 plus graphs (One complete set for each lab partner to be turned in; stapled together clearest report first.)


**DUE Next Lab Period** 

- (GQ) On-line *Density & Buoyancy Guiding Questions* (individually done)

  <a href="https://phet.colorado.edu/sims/density-and-buoyancy/density\_en.html">https://phet.colorado.edu/sims/density-and-buoyancy/density\_en.html</a>
- DUE Next Lab Period


### (GQ) On-line Density & Buoyancy Guiding DUE Next Lab Period


http://chemconnections.org/general/chem108/Density%20& %20Bouyancy.html



(GQ) On-line Density & Buoyancy Guiding Questions
 DUE Next Lab Period

https://phet.colorado.edu/sims/density-and-buoyancy/density\_en.html





Come to Dr. R. with partners having names on both lab forms.

Dr. R will provide unknown liquid, caliper & ruler. Metal cylinders are to be shared.

Read instructions carefully & discuss with partner before beginning.

Return empty liquid vials & metal cylinders when finished. Liquid is salt solution that can be poured down drain.